


PRODUCT DATA SHEET

Non-Functionalized Standard Gold Nanoparticles, >99% Endotoxin Free



Citrate Capped

Catalog Numbers: GEF-5-XX, GEF-10-XX, GEF-15-XX, GEF-20-XX, GEF-30-XX, GEF-40-XX, GEF-50-XX, GEF-60-XX, GEF-70-XX, GEF-80-XX, GEF-90-XX, GEF-100-XX, GEF-150-XX, GEF-200-XX, GEF-250-XX, GEF-300-XX, GEF-400-XX

Description

Cytodiagnosics Non-Functionalized Endotoxin Free Gold Nanoparticles have been optimized for cellular studies and high protein binding efficiency.

Our gold nanoparticles are available in 17 different sizes ranging from 5-400nm, are more than 95% spherical and have uniform size distribution (CV <12%).

For custom sizes, formulations or bulk quantities please contact our customer service department.

Features

- Superior size distribution compared to the leading competitor; available from 5nm to 400nm.
- Exceptional protein binding characteristics
- Ideal for sensitive applications requiring a minimal amount of contaminating reactants remaining from the manufacturing process. (Purity >99%)
- Endotoxin-free as tested by the Limulus amoebocyte lysate (LAL) test for endotoxins and determined to be <0.5 EU/mL, whilst also being free from substances that interfere with the endotoxin testing reaction.
- Packaged with dust-free protection to maintain the highest level of sterility, also making them suitable for use in clean rooms
- Also available with multiple surface functionalities to suit all your needs. (Please see Cytodiagnosics PEGylated and Functionalized Gold Nanoparticles)

Applications

- Excellent for development of peptide and protein gold conjugates for use in applications such as blotting, lateral flow assays, microscopy and transmission electron microscopy (TEM).
- Ideal for cellular and cellular toxicity studies, immunological studies, experimental models, and aseptic assays

Canada, Europe, Asia, Pacific and Africa

919 Fraser Drive, Unit 11, Burlington, ON Canada L7L 4X8

Tel: 866-344-3954 Fax: 289-204-9100

www.cytodiagnosics.com

Characteristics




Core diameter: 5-400nm (Coefficient of Variance < 12%)

Polydispersity Index (PDI): < 0.20

Concentration: OD=1 (~ 0.05 mg/ml)

Absorbance (λ_{max}): 510-570nm

Supplied in 0.1mM Phosphate-Buffered Saline (0.01X PBS)



Protein Conjugation

A recommended starting protocol for conjugation of proteins to Cytodiagnosics gold nanoparticles can be found online.

United States, Mexico, South and Central America

5867 South Garnett Road, Tulsa, Oklahoma 74146 USA

Tel: 866-344-3954 Fax: 289-204-9100

www.cytodiagnosics-us.com



Diameter (nm)	Peak SPR Wavelength (nm)	NPS/ml	Wt. Conc. (mg/ml)	Molar Ext (M ⁻¹ cm ⁻¹)	Size Dispersity (+/-nm)	Particle Volume (nm ³)	Surface Area (nm ²)	Surface/Volume Ratio	Particle Mass (g)	Molar Mass (g/mol)	Molar Conc.
5	515-520	5.47E+13	6.94E-02	1.10E+07	<15%	6.54E+01	7.85E+01	1.2	1.27E-18	7.64E+05	9.08E-08
10	515-520	5.98E+12	6.07E-02	1.01E+08	<15%	5.24E+02	3.14E+02	0.6	1.02E-17	6.11E+06	9.93E-09
15	520	1.64E+12	5.61E-02	3.67E+08	<12%	1.77E+03	7.07E+02	0.4	3.43E-17	2.06E+07	2.72E-09
20	524	6.54E+11	5.31E-02	9.21E+08	<12%	4.19E+03	1.26E+03	0.3	8.12E-17	4.89E+07	1.09E-09
30	526	1.79E+11	4.91E-02	3.36E+09	<12%	1.41E+04	2.83E+03	0.2	2.74E-16	1.65E+08	2.98E-10
40	530	7.15E+10	4.65E-02	8.42E+09	<12%	3.35E+04	5.03E+03	0.15	6.50E-16	3.91E+08	1.19E-10
50	535	3.51E+10	4.45E-02	1.72E+10	<10%	6.54E+04	7.85E+03	0.12	1.27E-15	7.64E+08	5.83E-11
60	540	1.96E+10	4.30E-02	3.07E+10	<10%	1.13E+05	1.13E+04	0.1	2.19E-15	1.32E+09	3.25E-11
70	548	1.20E+10	4.17E-02	5.03E+10	<10%	1.80E+05	1.54E+04	0.086	3.48E-15	2.10E+09	1.99E-11
80	553	7.82E+09	4.06E-02	7.70E+10	<10%	2.68E+05	2.01E+04	0.075	5.20E-15	3.13E+09	1.30E-11
90	564	5.37E+09	3.97E-02	1.12E+11	<8%	3.82E+05	2.54E+04	0.067	7.40E-15	4.46E+09	8.92E-12
100	572	3.84E+09	3.89E-02	1.57E+11	<8%	5.24E+05	3.14E+04	0.06	1.02E-14	6.11E+09	6.37E-12
150		3.60E+09	1.20E-01		<8%	1.77E+06	7.07E+04	0.04	3.43E-14	2.06E+10	5.98E-12
200		1.91E+09	1.55E-01		<8%	4.19E+06	1.26E+05	0.03	8.12E-14	4.89E+10	3.17E-12
250		7.08E+08	1.13E-01		<8%	8.18E+06	1.96E+05	0.024	1.59E-13	9.55E+10	1.18E-12
300		4.50E+08	1.24E-01		<8%	1.41E+07	2.83E+05	0.02	2.74E-13	1.65E+11	7.48E-13
400		1.88E+08	1.22E-01		<8%	3.35E+07	5.03E+05	0.015	6.50E-13	3.91E+11	3.12E-13

Storage

This product should be stored at 4°C. **DO NOT FREEZE.** If stored unopened and as specified, Cytodiagnosics gold nanoparticles are stable for at least 6 months.

Handling

When stored for a long period of time gold nanoparticles may sediment at the bottom of the flask, which is especially true for larger particle sizes. Prior to use, re-suspend the sedimented particles by swirling until a homogenous solution is obtained.

To maintain optimal performance, and stability of the colloidal gold, care should be taken to use clean storage containers if using other than supplied with the product.

Washing Gold Nanoparticles

Although it is not generally necessary to wash the gold nanoparticles prior to use, some applications may require additional washing procedures. The easiest way to remove possible contaminants in the nanoparticle solution is by centrifugation. Centrifugation force is dependent on size of

the gold nanoparticles and should be adjusted according to Table I for optimal performance.

Note I: Since non-functionalized gold nanoparticles are sensitive to salt containing buffers, re-suspension should always be performed in ultra-pure water to prevent irreversible aggregation. Irreversible aggregation is characterized by a clear to bluish solution upon the addition of salt.

Note II: Please note that centrifugation can induce aggregation. To prevent aggregation, it may be necessary to add Tween 20 at a concentration of 0.025% w/v.

Canada, Europe, Asia, Pacific and Africa

919 Fraser Drive, Unit 11, Burlington, ON Canada L7L 4X8
Tel: 866-344-3954 Fax: 289-204-9100
www.cytodiagnosics.com

United States, Mexico, South and Central America

5867 South Garnett Road, Tulsa, Oklahoma 74146 USA
Tel: 866-344-3954 Fax: 289-204-9100
www.cytodiagnosics-us.com



Procedure

1. Place 1ml aliquot of colloidal gold in a 1.5ml micro centrifuge tube.
2. Centrifuge the gold nanoparticles for 30 minutes using the appropriate G force determined by referencing Table I.
3. Remove the supernatant and re-suspend in an appropriate volume of ultra-pure water.
4. Vortex to re-disperse the particles.

Product Safety and Handling

This product is for R&D use only, not for drug, household, or other uses. Please review the material safety datasheet (MSDS) available online for proper safety and handling procedures.

Ordering Information

For ordering call 866-344-3954 or visit us online.

Table I. Appropriate G forces for centrifugation of gold nanoparticles. Note that recommended conditions are for a volume of 1ml and centrifugation using a microcentrifuge,

Size (nm)	Speed (g)	Time (min)
5	100,000	30
10	17,000	60 (~50% recovery)
15	17,000	30
20	6,500	30
30	4,500	30
40	2,500	30
50	2,000	30
60	1,125	30
80	600	30
100	400	30
150	180	30
200	100	30

Canada, Europe, Asia, Pacific and Africa

919 Fraser Drive, Unit 11, Burlington, ON Canada L7L 4X8
Tel: 866-344-3954 Fax: 289-204-9100
www.cytodiagnosics.com

United States, Mexico, South and Central America

5867 South Garnett Road, Tulsa, Oklahoma 74146 USA
Tel: 866-344-3954 Fax: 289-204-9100
www.cytodiagnosics-us.com



Catalogue Number	Description	Sizes
GEF-5-XX*	5nm Endotoxin Free Gold Nanoparticles	20ml, 100ml, 500ml or 1000ml
GEF-10-XX	10nm Endotoxin Free Gold Nanoparticles	20ml, 100ml, 500ml or 1000ml
GEF-15-XX	15nm Endotoxin Free Gold Nanoparticles	20ml, 100ml, 500ml or 1000ml
GEF-20-XX	20nm Endotoxin Free Gold Nanoparticles	20ml, 100ml, 500ml or 1000ml
GEF-30-XX	30nm Endotoxin Free Gold Nanoparticles	20ml, 100ml, 500ml or 1000ml
GEF-40-XX	40nm Endotoxin Free Gold Nanoparticles	20ml, 100ml, 500ml or 1000ml
GEF-50-XX	50nm Endotoxin Free Gold Nanoparticles	20ml, 100ml, 500ml or 1000ml
GEF-60-XX	60nm Endotoxin Free Gold Nanoparticles	20ml, 100ml, 500ml or 1000ml
GEF-70-XX	70nm Endotoxin Free Gold Nanoparticles	20ml, 100ml, 500ml or 1000ml
GEF-80-XX	80nm Endotoxin Free Gold Nanoparticles	20ml, 100ml, 500ml or 1000ml
GEF-90-XX	90nm Endotoxin Free Gold Nanoparticles	20ml, 100ml, 500ml or 1000ml
GEF-100-XX	100nm Endotoxin Free Gold Nanoparticles	20ml, 100ml, 500ml or 1000ml
GEF-150-XX	150nm Endotoxin Free Gold Nanoparticles	20ml, 100ml, or 500ml
GEF-200-XX	200nm Endotoxin Free Gold Nanoparticles	20ml, 100ml, or 500ml
GEF-250-XX	250nm Endotoxin Free Gold Nanoparticles	20ml, 100ml, or 500ml
GEF-300-XX	300nm Endotoxin Free Gold Nanoparticles	20ml, 100ml, or 500ml
GEF-400-XX	400nm Endotoxin Free Gold Nanoparticles	20ml, 100ml, or 500ml

*Indicates quantity, e.g., GEF-5-20 for 20ml of 5nm Endotoxin Free Gold Nanoparticles

For custom sizes (up to 400nm available), bulk quantities, and custom gold nanoparticle surface chemistry, please contact our customer service department.

Canada, Europe, Asia, Pacific and Africa

919 Fraser Drive, Unit 11, Burlington, ON Canada L7L 4X8
Tel: 866-344-3954 Fax: 289-204-9100
www.cytodiagnosics.com

United States, Mexico, South and Central America

5867 South Garnett Road, Tulsa, Oklahoma 74146 USA
Tel: 866-344-3954 Fax: 289-204-9100
www.cytodiagnosics-us.com