

miriad®

RVF Toolkit



Features

- Instant assay results
- Flexible sample types and volumes
- No specialized equipment required
- Simultaneous detection of multiple biomarkers on a single device
- Room-temperature storage
- Versatility for multiple applications in laboratories or in the field
- Professional consultation and licensing opportunities with MedMira

Rapid assays reinvented.

MedMira's Miriad® RVF Toolkit is an off-the-shelf package of components and protocols. The toolkit provides the building blocks needed to create and optimize assays whose performance can rival much more complex laboratory-based systems.



These assays can be used in developing rapid tests to detect an analyte of interest, to monitor bioreactors or be applied to other processes in a commercial lab. The toolkit can also provide a simple migration route for existing ELISA and Western blot assays to become rapid assays that can be used in-house or commercialized.

Quickly create rapid assays in your lab.

Use MedMira's RVF Toolkit, plus materials and supplies in your lab, to create rapid analytical tools in three steps.

1. Spot your antigen
2. Let it dry
3. Run your test

Potential Research Applications and Opportunities

Infectious Disease Biomarkers

- Demonstrate the functionality of an antigen
- Produce a sensitive and specific test for commercial use
- Use for human or veterinary purposes
- Monitor early stage experiments

Vaccine Research

- Create alpha prototypes employing antigens used in vaccine development
- Determine the ability of a vaccine to illicit immune responses
- Monitor methods used to produce the antigen

Industrial Processes


- Develop rapid prototypes
- Monitor bioreactor production of various antibodies
- Use in a variety of applications, including food safety and quality control

Contamination Identification

- Seek contaminants in cell culture, DNA or RNA
- Create single or multiplexed rapid test to diagnose issue in real-time

ELISA and Western Blot Alternatives

Migrate labor-intensive and long-running ELISA and Western blot assays to the Rapid Vertical Flow platform quickly and cost-effectively



Rapid Vertical Flow Technology®

At the core of the Miriad® toolkit is MedMira's patented Rapid Vertical Flow technology®, offering speed, accuracy and multiplexing. This technology enables targets of interest in a wide variety of sample types to be captured and visualized on a unique membrane. MedMira's InstantGold™ Cap simplifies the visualization procedure, as well as the handling and storage of the detection agent. Rapid Vertical Flow technology® is a highly scalable solution that is transforming the development and commercialization of rapid testing applications across a broad range of life sciences and biotech sectors.

Technical Support and Consultation

MedMira offers technical assistance, consulting services and licensing agreements. These services can help optimize assays, as well as support the validation and commercialization of discoveries.

We are experienced with a variety of applications such as:

- Assay development
- Prototype development
- Manufacturing of rapid based immunoassays
- Preparation, documentation, and maintenance of design plans and design history files
- Laboratory studies
- Analytical sensitivity and specificity studies
- Expert regulatory guidance on all design and manufacturing issues
- Manufacturing product post-approval

About MedMira

MedMira's patented Rapid Vertical Flow technology delivers fast, accurate results in a range of applications. Our products help care providers and their patients know more, instantly. MedMira rapid tests are manufactured in a cGMP facility with a quality management system registered to ISO 13485:2016. Learn more at medmira.com

Ordering Information

Cat. No.: RVF-05-50

Description: Miriad®
RVF Toolkit

Quantity Box: 50



MedMira Laboratories Inc.
155 Chain Lake Drive, Suite 1
Halifax, NS B3S 1B3 CANADA
T. | 902.450.1588
E. support@medmira.com



For Canada, Europe, Africa,
Asia/Pacific
Distributor:
Cytodiagnosics Inc.
919 Fraser Dr., Unit 11
Burlington, ON, L7L 4X8
Canada
orders@cytodiagnosics.com
Phone: 866-344-3954

For USA, Mexico, Central
America, South America
Distributor:
Cyto US Inc.
5867 South Garnett Road
Tulsa, Oklahoma, 74146
United States
orders-us@cytodiagnosics.com
Phone: 866-344-3954

